

Northern Blot (sRNA) -- LICOR

General Notes:

Practice RNase-Free technique:

- Use RNA-only pipet tips/tubes
- Use RNaseZap to clean counters and pipetmen
- Filter/sterilize solutions
- Use DEPC-H₂O or MQ-H₂O
- Keep RNA cold (4 °C or lower)

Needed Material:

- Mini-PROTEAN Tetra Cell
- Mini Trans-Blot Cell Assembly for Transfer
- Loading buffer (LBI (Ambion)- Glycerol, bromophenol blue, xylene cyanol, formamide), stored at -20 °C

- 10 X TBE Stock (Use 1 X TBE dilution – F/S):

- 108 g Tris Base
- 55 g Boric Acid
- 40 mL 0.5 M EDTA (disodium salt), pH 8.0
- In 960 mL MQ H₂O

- 20 X SSC (1 X and 6 X dilutions – F/S):

- 1X: 1 mL 20X SSC
- 19 mL MQ H₂O

- 6X: 6 mL 20X SSC
- 14 mL MQ H₂O

-ULTRAhyb-Oligo (Ambion), stored at 4 °C, needs to be heated to go into solution prior to use

-Odyssey Blocking Buffer (LiCor)

-Wash Buffers (does not need to be filter sterilized)

- Low: 50 mL 20X SSC
- 2.5 mL 20% SDS
- Bring to 500 mLs with MQ H₂O

- High: 2.5 mL 20X SSC
- 2.5 mL 20% SDS
- Bring to 500 mLs with MQ H₂O

Protocol (for detecting ~120 nt sRNA):**1. Prepare a 10% TBE-Urea gel:**

To a beaker **add**:

- a. 3.3 mL 30% acrylamide (protogel)
- b. 1.0 mL 10 X TBE
- c. 2.136 mL water (MQ)
- d. 4.2 g urea

Cover with Parafilm and Incubate at 37°C or in hot bath (65°C) until urea dissolves. Then add

- e. 330 µL 1.6% APS (Needs to be less than 1 month old)
- f. 5 µL TEMED

Swirl to facilitate polymerization and then pour into gel casting frame.

Allow gel to solidify in a Mini-PROTEAN Tetra cell casting frame/stand for 45'. Gel can be stored overnight, wet, at 4 °C.

2. Prepare RNA Samples:

- a. **Prepare** RNA samples with equal volume of LBII, (try to load the same volume in each lane; adjust volume of RNA with DEPC-H₂O).
- b. **Heat** for 5' at 85 °C, **Pico-centrifuge**.
- c. Keep RNA cold (4 °C or lower).

3. Pre-Run and Run Gel:

- a. **Pre-Run gel** in 1 X TBE at 250 V for 20' in a Mini-PROTEAN Tetra cell tank.
- b. **Pipette** out gel debris from each lane before loading RNA samples.
- c. **Run gel** in 1 X TBE at 200 V for 45-60', (ensure dark blue loading dye runs off gel).

****If you are planning to use the Trans Blot Turbo, begin preparation of Riboprobe (see riboprobe protocol)****

4. Transfer Gel to Membrane:**a. Traditional Method**

- i. **Equilibrate** sponges, filter paper and membrane in 1x TBE.
- ii. **Prepare** gel sandwich (black, sponge, filter paper, 10% TBE-Urea gel, a 60 cm² membrane, filter paper, sponge, white).
- iii. **Transfer** (black to black, clear to red) in Mini Trans-Blot Cell in 1 X TBE at 200 mA for 60'. Add a stir bar and ice block to cell and place on a magnetic stirrer.

****Begin preparation of Riboprobe (see riboprobe protocol)****

b. Trans Blot Turbo

- i. **Equilibrate** extra-thick filter paper and membrane in 1x TBE.
- ii. **Prepare** gel sandwich in the Trans-Blot Turbo cassette B (cassette bottom, extra-thick filter paper, membrane, TBE-Urea gel, extra-thick filter paper, cassette top)
- iii. **Transfer** using turbo command on main menu (third button from the left on the bottom row.) Select the number/size gel/s that you will be transferring with this program. This will take you to a screen showing the current conditions in the cassettes. Push the B:RUN button to start the transfer. Transfer will be finished in 7 minutes.

5. **Wash Membrane and cross-link:**

Handle membrane gently, avoid touching with fingers:

- a. **Wash membrane** in 20 mL of 6 X SSC for 2' with rocking.
- b. UV **Cross link: "Autocross link"**.
- c. **Wash membrane** in 20 mL of 1 X SSC for 1' with rocking.

6. **Prehybridization and Hybridization:**

- a. **Pre-hybridize** ($\geq 30'$) standard membrane size in 6 mL of ULTRAhyb-Oligo (Ambion) at 60-65 °C.
- b. Add biotinylated riboprobe specific for RNA (10 μ L of a 20 μ L in vitro transcription reaction is PLENTY!) and DNA probe for 5S (2.5-5 μ L, ~ 250 ng, 800-5S).
- c. Allow to **hybridize** overnight at 60-65 °C.

Following day, wash membrane using Li-Cor detection method. See Li-Cor Northern protocol.